2 × PCR Solution

PrimeSTAR™ HS (Premix)

Code No. R040A

Shipping at – 20°C

Stored at – 20°C

Size: 500 µl × 5
(for 100 PCR reactions)

Storage: Excessively repeated freeze-thaw cycles may decrease the enzyme activity.

Lot No.

Expiry Date:

Description: PrimeSTAR™ HS (Premix) is an optimized mixture composed of PrimeSTAR™ HS DNA Polymerase, which is a high-fidelity DNA polymerase developed by Takara Bio Inc, reaction buffer and dNTP mixture as 2-fold concentration. As this product offers quick preparation of reaction mixture and reduction of contamination risk, it is also useful for high-throughput application.

PrimeSTAR™ HS DNA Polymerase has a matchless proof reading activity due to very strong 3' -5' exonuclease activity, and besides its amplification efficiency is higher than that of Taq DNA Polymerase.

Content:

PrimeSTAR™ HS DNA Polymerase * : 1.25 units/ 25 µl
dNTP Mixture : 2 × conc. ; ea. 0.4 mM
PrimeSTAR™ Buffer : 2 × conc. ; including 2 mM Mg²⁺

* Specification of PrimeSTAR™ HS DNA Polymerase (Cat.#R010A)

Unit definition: One unit is the amount of the enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

100 mM	Tris-HCl (pH8.3 at 37°C)
10 mM	KCl
2 mM	MgCl ₂
6 mM	(NH ₄) ₂ SO ₄
0.1%	TritonX-100
200 µM	each dATP, dGTP, dCTP
100 µM	[³ H] TTP
0.001%	BSA
0.4 mg/ml	activated salmon sperm DNA

Purity: Nicking, endonuclease and exonuclease activity were not detected after incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ-Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

Applications: DNA amplification by Polymerase Chain Reaction (PCR).

PCR product: A significant percentage of PCR product obtained using PrimeSTAR™ HS DNA polymerase will possess blunt-ends. Thus, obtained PCR products can be directly cloned into blunt-end vectors. (If necessary, phosphorylate PCR products before cloning.)

PCR test: Good enzyme performance was confirmed by robust single fragment DNA PCR amplification using both λ DNA (amplified fragments: 8, 10, 12, 15 kb) and Human genomic DNA (amplified fragment: 0.5, 1, 2, 4, 6, 8 kb) as templates.

General reaction mixture for PCR (total 50 µl):

PrimeSTAR™ HS (Premix)	25 µl
Primer 1	0.2-0.3 µM (final conc.)
Primer 2	0.2-0.3 µM (final conc.)
Template	<200 ng
Sterilized distilled water	up to 50 µl

Recommended template amount

Human genomic DNA	5 ng-200 ng (<200 ng)
<i>E. coli</i> genomic DNA	100 pg-100 ng
cDNA library	1-200 ng
λ DNA	10 pg-10 ng
plasmid DNA	10 pg-1 ng

PCR conditions:

98°C	10 sec.	30 cycles
55°C	5 sec. or 15 sec.	
72°C	1 min./kb	
or		
98°C	10 sec.	30 cycles
68°C	1 min./kb	

Notice: Because this enzyme possesses an extremely high priming efficiency, use of a short annealing time is strongly recommended. Please refer to the product manual supplied with this product for detailed cycling recommendations.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as the patented 5' Nuclease Process claims in US Patents Nos. 5,210,015 and 5,487,972), no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

[L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287 and corresponding patents in other countries.

[M54] PrimeSTAR™ HS DNA Polymerase

This product is the subject of the pending U.S. patent application and its foreign counterparts.

U.S. Patent 5,436,149 for LA Technology is owned by TAKARA BIO INC.

Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not be used for products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at www.takara-bio.com.

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